

# Biotage


*Please see Amanda for training before your first use.*

## Sample Preparation


1. Concentrate your sample into a solid or an oil. (You can also dry load your sample onto silica and pack it into one of the SNAP columns.)
2. Run a TLC on your compound that needs purification and determine its R<sub>f</sub> value along with the R<sub>f</sub> values of all impurities.

## Running a Sample


1. Press **Lamp On**. (Note that the lamp has a 7.5 min warm-up time.)
2. Select **Solvents** tab and make sure the solvents that you want to use are filled. Check also to see how full the waste bottle is at this time.
3. Then select the **Method** tab. At the bottom of the screen select **From TLC...** Enter your pre-determined TLC conditions along with the R<sub>f</sub> values you calculated. Select the appropriate column size you are using at this time and press close.



4. At the top of the screen open the **Parameters** tab, select the solvents, cartridge type, and rack type to be used. Enter a sample name for the run.



- At the **Collection** tab, select the UV wavelengths and start threshold.



- At the **Gradient** tab, review the elution gradient and if needed, drag the blue line to modify the gradient. (Select the **equilibrate** button if elution of the column is desired.)
- Load the cartridge and racks onto the system. Ensure that the tubing is securely attached.
- Press **Run**.
- In the **Run Parameters** dialog box, select the appropriate rack positions.
- To start the run or equilibration, press **Play**.
- After equilibration is done. Select the **Load sample** button and inject your sample on to the column. Afterwards, select **Gradient** to run the sample.

## Processing the Data

- When the run is complete all the data will be stored under the **Results** tab.
- In the **Results** tab select the different fractions collected, the gradient used, and the time of the entire run.